
Snap 40

Portable Alcohol Meter for Distillates

Portable Density & Concentration Meter

- ☒ Check with distilled water and alcohol at different temperatures.
- ☒ Check of all labels.

Yours sincerely,

A handwritten signature in blue ink, appearing to be 'Kle', is written over a horizontal dashed line.

Inspector

Certificate of Compliance

Anton Paar

For the

Snap 40

Serial Number: 82372286

Dear customer,

This instrument was inspected, tested and accepted in accordance with our quality system and the applicable specifications.

The following tests have been successfully performed:

- ☒ Optical check of electronic boards, wiring system and surface finish.
- ☒ Functional test of the instrument, including display and keyboard, pipette-style pump, automatic "power off" and acoustic signal
- ☒ Test of interfaces
- ☒ Measurements: Adjustment with air at 20 °C and distilled water at 8 °C, 20 °C, 40 °C

Contents

1 About the Instruction Manual	7
2 Safety Instructions	8
3 Measuring Principle	10
4 Snap 40 – an Overview	11
5 Checking the Supplied Parts	12
6 Description of the Instrument	14
6.1 Front View	14
6.2 Rear View	16
6.3 Display	17
7 Operating the Snap 40	18
7.1 Connecting the Filling Tube	18
7.2 Attaching the Syringe Adapter	19
7.3 Switching the Instrument On	19
7.4 Switching the Instrument Off	19
8 Defining the Basic Settings	20
8.1 Language Settings	20
8.2 Sound Settings	20
8.3 Energy Saving Mode	20
8.4 Backlight	20
8.5 Contrast Settings	21
9 Performing a Measurement	22
9.1 Selecting the Measuring Unit	23
9.1.1 Measuring Units	23
9.2 Filling Sample	23
9.2.1 Filling via the Filling Tube	24
9.2.2 Filling via the Plastic Syringe	25
9.3 Performing a Measurement	25
9.3.1 Emptying the Measuring Cell	25

10 Cleaning and Storing the Instrument	26
10.1 Cleaning the Measuring Cell	26
10.2 Cleaning the Filling Pump	27
10.3 Cleaning the Pump Housing	27
10.4 Cleaning the Housing and Display	27
10.5 Storing the Instrument	27
11 Maintenance	28
11.1 Changing the Batteries	28
11.2 Carrying out a Firmware Update	30
11.3 Device Information	31
12 Readjustment	32
12.1 Check Measurement	32
12.2 Water Adjustment	32
12.3 Restoring to Factory Adjustment	33
13 Error Messages and Warnings	34
14 Technical Data	35
15 Wetted Parts	36
Appendix A: Density of Water (0 °C to 40 °C)	37
Appendix B: Measuring Unit Details	38
Appendix C: CE Declarations of Conformity	39
Appendix D: Firmware Versions	40
Appendix E: Menu Tree	40

Snap 40

Portable Alcohol Meter for Distillates

From Firmware Version V1.00

Instruction Manual

While every precaution has been taken in the preparation of this document, Anton Paar GmbH assumes no responsibility for technical or printing errors or omissions. Nor is any liability assumed for damages resulting from the use of the information contained in this instruction manual. Anton Paar GmbH does not make a commitment to update the information in this manual. Specifications are subject to change without notice.

All rights reserved (including translation). No part of this document may be translated, reproduced or distributed in any form (print, photocopy, microfilm or any other process) without the prior written permission of Anton Paar GmbH. Trade marks may be used in this instruction manual without being marked as such. These are the property of their respective owners and are legally protected.

Published by Anton Paar GmbH. Printed in Austria. Copyright © 2014 Anton Paar GmbH, Graz, Austria

Contact: Anton Paar GmbH
Anton-Paar-Str. 20
A-8054 Graz / Austria - Europe
Tel: +43 316 257-0
Fax: +43 316 257-257
E-mail: info@anton-paar.com
Web: www.anton-paar.com

Date: May 06, 2014
Document number: C96IB005ML-A

1 About the Instruction Manual

This instruction manual informs you about the installation and the safe handling and use of the product. Pay special attention to the safety instructions and warnings in the manual and on the product.

The instruction manual is a part of the product. Keep this instruction manual for the complete working life of the product and make sure it is easily accessible to all people involved with the product.

Conventions for safety messages

The following conventions for safety messages are used in this instruction manual:

DANGER

Danger indicates a hazardous situation which, if not avoided, will result in death or serious injury.

WARNING

Warning indicates a hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION

Caution indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.

NOTICE

Notice indicates a situation which, if not avoided, could result in damage to property.

TIP *Tip gives extra information about the situation at hand.*

Typographical conventions

The following typographical conventions are used in this instruction manual:

Convention	Description
<key>	The names of keys and buttons are written inside angle brackets.
"Menu Level 1 > Menu Level 2"	Menu paths are written in bold, inside straight quotation marks. The menu levels are connected using a closing angle bracket.

2 Safety Instructions

- Read this instruction manual before using the portable alcohol meter Snap 40.
- Follow all hints and instructions contained in this instruction manual to ensure the correct use and safe functioning of Snap 40.

Liability

- This instruction manual does not claim to address all safety issues associated with the use of the instrument and samples. It is your responsibility to establish health and safety practices and determine the applicability of regulatory limitations.
- Anton Paar GmbH only warrants the proper functioning of Snap 40 if no adjustments have been made to the mechanics, electronics, and firmware.
- Only use Snap 40 for the purpose described in this instruction manual. Anton Paar GmbH is not liable for damages caused by incorrect use of Snap 40.

Installation and use

- Snap 40 is **not** explosion-proof instruments and therefore must not be operated in areas with risk of explosion.
- Never remove the pump lock and battery cover in hazardous areas. Only exchange the batteries outside of hazardous areas.
- The installation procedure shall only be carried out by authorized personnel who are familiar with the installation instructions.
- Do not use any accessories or wearing parts other than those supplied or approved by Anton Paar GmbH.

- Make sure all operators are trained to use the instrument safely and correctly before starting any applicable operations.
- In case of damage or malfunction, do not continue operating Snap 40. Do not operate the instrument under conditions which could result in damage to goods and/or injuries and loss of life.
- Check Snap 40 for chemical resistance to the samples and cleaning liquids.
- Do not expose the instrument to temperatures below 0 °C (32 °F) if water is contained in the measuring cell or pump (freezing water will cause rupture of the measuring cell).

Maintenance and service

- The results delivered by Snap 40 not only depend on the correct functioning of the instrument, but also on various other factors. We therefore recommend you have the results checked (e.g. plausibility tested) by skilled personnel before consequential actions are taken based on the results.
- Service and repair procedures may only be carried out by authorized personnel or by Anton Paar GmbH.
- If your instruments needs repair, contact your local Anton Paar representative.

Disposal

- Concerning the disposal of Snap 40 observe the legal requirements in your country.

Precautions for highly flammable samples and cleaning liquids

- Observe and adhere to your national safety regulations for handling the measured samples (e.g. use of safety goggles, gloves, respiratory protection etc.).
- Only store the minimum required amount of sample, cleaning liquids and other inflammable materials near Snap 40.
- Do not leave sample/rinsing containers uncovered. Clean all spillages immediately.
- Make sure that the setup location is sufficiently ventilated. The environment of Snap 40 must be kept free of flammable gases and vapors.
- Supply a fire extinguisher.

3 Measuring Principle

Snap 40 determines the alcohol concentration of a distillate based on its density.

Definition of density

The density ρ of a sample is defined as its mass divided by its volume:

$$\rho = \frac{m}{V}$$

Density is a temperature-dependent measuring unit.

The oscillating U-tube principle

The sample is introduced into a U-shaped borosilicate glass tube that is being excited to vibrate at its characteristic frequency electronically. The characteristic frequency changes depending on the density of the sample. Through determination of the characteristic frequency the density of the sample can be calculated. Due to the temperature dependency of the density value, the temperature of the sample has to be determined precisely.

Further the temperature of the sample should be close to the ambient temperature and instrument's temperature to keep the temperature influence on the measured result to a minimum. Density as well as temperature are measured simultaneously and are the basis for the calculation of the % alcohol, compensated to a reference temperature.

Concentration measurement

In binary mixtures, the density of the mixture is a function of its composition. Thus, by using density/concentration tables, the density value of a binary mixture can be used to calculate its composition.

This is also possible with so-called quasi binary mixtures. These are mixtures containing two major components and some additional ones which are present in very small concentrations compared to the two main components. Many decarbonated soft drinks, for example, can be considered to be quasi binary mixtures of sugar in water because the concentration of flavors and acids are very small compared to sugar and water. Thus, the sugar concentration can be measured with a density meter.

The same statement is true for measurement of the alcohol concentration in distillates that can be considered to be binary mixtures of ethanol and water.

In contrast, all kind of liqueurs have to be considered as ternary solutions of water, ethanol and sugar. They are typically produced by mixing distillates with fruit juices, herbs, dairy products etc. and adding sugar. Therefore the determination of the true alcohol concentration based on density measurement is not possible for liqueurs.

4 Snap 40 – an Overview

Snap 40 is a portable alcohol meter for determination of the alcohol concentration in distillates that can be considered to be binary mixtures of ethanol and water. The concentration determination is based on density measurement according to the oscillating U-tube principle.

Quick on-site measurements

The robust instrument measures samples directly out of the sample container – without the necessity to pour the sample into a measuring cylinder. After pressing down the lever of the pump, the tube is immersed directly into the storage vessel to suck up 2 mL of sample. Within seconds the instrument displays the alcohol concentration with an accuracy of 0.2 % v/v.

One instrument for distillates of all strengths

Snap 40 covers a measuring range from 0 % v/v to 100 % v/v (0 to 200 °Proof) and therefore measures distillates of all strengths – be it pure distillates, samples during dilution of pure distillate to drinking strength or a check of the final product.

Automatic temperature compensation

Snap 40 measures sample temperatures from 5 to 30 °C and automatically compensates the measured result to the required reference temperature (20 °C for % v/v alcohol, 60 °F for °Proof). The final result is clearly visible on the large display.

Easy operation

The measuring procedure is simple and the user interface is so intuitive that you will be familiar with the instrument from the start. On top of this, Snap 40 is easy to clean, robust and ready for a long working life.

5 Checking the Supplied Parts

Snap 40 was tested and packed carefully before shipment. However, damage may occur during transport.

1. Keep the packaging material (box, foam piece) for possible returns and further questions from the transport company or insurance company.
2. Check the delivery for completeness by comparing the supplied parts to the list in Table 5.1.
3. If a part is missing, contact your Anton Paar representative.
4. If a part is damaged, contact the transport company and your Anton Paar representative.

Table 5.1: Supplied parts

Symbol	Pcs.	Article Description	Mat. No.
	1	Snap 40 portable alcohol meter	135166
	1	Instruction Manual English/German	146924
	1	Filling tube (standard 180 mm)	68527
	1	Luer adapter 1/4" UNF (for syringe filling)	64792
	10	Plastic syringes 2 mL	58802
	1	Allen wrench 2.5 mm DIN 911	58263

Table 5.2: Optional parts

Article Description	Mat. No.
Carrying case	109048
Filling tube PTFE, length: 600 mm	78503
IrDA USB adapter LCS-8141	88085
Set wristband for portable instruments	92416
Rubber housing DMA 35/Snap 40	105540

6 Description of the Instrument

6.1 Front View

Fig. 6 - 1 Front view of Snap 40

- 1 ... Built-in pump
- 2 ... Screw plug
- 3 ... Measuring cell
- 4 ... Softkeys
- 5 ... Operating keys
- 6 ... Graphical, monochrome LC display


Fig. 6 - 2 Top view of Snap 40

- 1 ... Built-in pump
- 2 ... Fixing screw of the pump lock
- 3 ... Pump lock
- 4 ... Infrared interface (IrDA)

Keys at the front side

	For switching the instrument on and off.
	For starting the water check.
	Softkeys for selecting menu items and for navigation. The function of the right softkey can be configured.
	Arrow keys for navigation within the menu and for the entry of characters.

TIP For faster up and down navigation, keep the arrow keys pressed.

6.2 Rear View


Fig. 6 - 3 Rear view of Snap 40

- 1 ... Type plate with serial number
- 2 ... Mark for the correct battery insertion

6.3 Display

Fig. 6 - 4 Display of Snap 40

- | | |
|--|--------------------------------|
| 1 ... Information header | 5 ... Measuring value |
| 2 ... Status of measurement | 6 ... Sample temperature °C/°F |
| 3 ... Symbol for battery charge status | 7 ... Softkey configurations |
| 4 ... Warning sign | |

Symbols on the display

	The infrared data transfer symbol is shown on the information header during a firmware update via the IrDA interface.
	The battery charge status symbol informs you about the battery charge status of your instrument. For further information on different charge status symbols and their meanings, see Chapter 11.1.
	The warning sign indicates that a measured value is outside range specifications. The warning sign is always shown together with an information (press the right softkey) stating the warning type (for details about warning types, see Chapter 13).

7 Operating the Snap 40

NOTICE

Standard tolerances in dimensions of AA batteries might result in some rattling noise coming from the battery compartment when moving the instrument. This has no negative influence on the operation of the instrument.

NOTICE

Protect the instrument against any collision or strong shock as this may lead to a malfunction of the instrument and breakage of the measuring cell.

7.1 Connecting the Filling Tube

- Screw in the filling tube by hand until you feel some resistance against turning. Do not use any tools for screwing in the filling tube.

Fig. 7 - 1 Connecting the filling tube

7.2 Attaching the Syringe Adapter

1. Remove the screw plug at the side of the instrument.
2. Screw in the Luer adapter (see Chapter 5) by hand until you feel some resistance against turning. Do not use any tools for screwing in the Luer adapter.

7.3 Switching the Instrument On

- To switch the instrument on, press the  key for approximately two seconds until the instrument beeps and shows the welcome screen.

After showing the welcome screen, Snap 40 immediately switches to the main screen (measuring mode). Now you can start with the measuring process.

7.4 Switching the Instrument Off

- To switch the instrument off, press the  key until "Power Off" is displayed.

If the instrument is in the energy saving mode (see Chapter 8.3), Snap 40 switches automatically off after 3, 5 or 10 minutes inactivity (depending on the setting).

8 Defining the Basic Settings

8.1 Language Settings

You can choose between four different menu languages:

1. Press <Menu> and select **"Setup > Language"**.
2. Select the preferred language (English, German, French or Spanish).
3. Press <OK>.

8.2 Sound Settings

When sound is enabled, Snap 40 beeps when it is switched on and off and when a key is pressed.

When sound is disabled, the instrument only beeps when it is switched on and off.

1. Press <Menu> and select **"Setup > Beep"**.
2. Select the preferred setting (Beep on/off) and press <Save>.

8.3 Energy Saving Mode

Enabling the energy saving mode lowers the power consumption of the instrument and thus extends battery life. When the energy saving mode is enabled, the instrument will automatically switch itself off after 3, 5 or 10 minutes of inactivity.

When the energy saving mode is disabled, the instrument will remain switched on until it is turned off using the  key.

1. Press <Menu> and select **"Setup > Energy Saving"**.
2. Select the preferred setting (Off, 3, 5 or 10 minutes) and press <Save>.

8.4 Backlight

The LC display and the oscillator are equipped with a backlight. The backlight ensures that even under bad lighting conditions the measuring results and menu options are easily readable. The backlight of the oscillator enables you to observe the filling process.

You can choose between "Auto", "Off" and "On":

Auto	The backlight automatically turns on when a key is pressed or when a new sample is filled into the measuring cell. To activate the backlight without accessing the menu, press an arrow key. The backlight switches off automatically after 30 seconds.
Off	The backlight is permanently off.
On	The backlight is permanently on.

TIP Having the backlight turned on permanently will shorten battery life.

1. Press <Menu> and select **"Setup > Backlight"**.
2. Select the preferred setting (Auto, Off, On) and press <Save>.

8.5 Contrast Settings

The display contrast can be freely adjusted to suit your preference.

1. Press <Menu> and select **"Setup > Display Contrast"**.
2. Select the preferred setting (-8 to +8) and press <Save>.

9 Performing a Measurement

General instructions for measuring

WARNING

Handling samples with temperatures of more than 70 °C bears the danger of heavy burns.

- Make sure you wear protective clothes or ensure alternative protection from burns when handling high temperature samples.

NOTICE

Before performing a measurement, make sure that the wetted parts are resistant to the sample and cleaning liquids (see Chapter 15).

NOTICE

Sample containing dissolved CO₂ will cause bubbles within the measuring cell with the effect of invalid measurement results. Make sure your sample is degassed carefully.

Degassing a sample can be done by:

- boiling the sample for several minutes,
 - stirring the sample for 5 to 15 minutes until no bubbling occurs any more, or
 - putting the sample for approximately 5 to 10 minutes into an ultrasonic bath.
-
- Make sure that the measuring cell is filled free from gas bubbles. Possible reasons for gas bubbles in the measuring cell are:
 - Gas bubbles in the sample
 - Leaky connection of the filling tube, the pump or the screw plug
 - Make sure that the measuring cell is fully filled.
 - Make sure that the sample temperature does not deviate too strongly from the ambient temperature.
- If the sample measured is outside the temperature specifications (5 to +30 °C), the display shows the warning sign until the sample has reached a measurable temperature.
- Make sure to carefully clean the instrument after each measurement series to avoid deposits in the measuring cell.
 - Make sure that suitable solvents for cleaning are available.

TIP Your filled sample is continuously measured and the measuring result is displayed.

9.1 Selecting the Measuring Unit

For measuring your samples you can choose between two predefined measuring units.

Both measuring units are derived from the density of the sample at the measured temperature.

1. Press <Menu> and select "Units".

The list of the predefined measuring units is displayed.

2. Select the preferred measuring unit and press <Save>.

9.1.1 Measuring Units

For further details, see Appendix B.

Meas. units	Options	Description
Alcohol	Alcohol % v/v @ 20 °C	Concentration of an ethanol/water mixture in % by volume at 20 °C.
	Alcohol US @ 60 °F ("Proof")	Degrees Proof at 60 °F.

9.2 Filling Sample

You can fill the measuring cell using the filling tube or the plastic syringe. When filling highly viscous samples, we recommend using the plastic syringe.

9.2.1 Filling via the Filling Tube

1. Press down the pump lever as far as it will go (see Fig. 9 - 2).
2. Submerge the filling tube in the sample.


Fig. 9 - 1 Pressing down the pump lever and submerging the filling tube

3. Slowly release the pump lever.

Fig. 9 - 2 Releasing the pump lever and reading the result

9.2.2 Filling via the Plastic Syringe

NOTICE

Do not remove the pump when filling with the plastic syringe. Otherwise, the system is leaky.

1. Fill the plastic syringe with the sample.
2. Fill the measuring cell through the Luer adapter using the plastic syringe.

Fig. 9-3 Filling sample using the plastic syringe

9.3 Performing a Measurement

1. Select the measuring unit (see Chapter 9.1).
2. Fill the measuring cell with the sample.
The measured value is displayed continuously. In the information header the measuring status is displayed. Wait until the status turns "Valid" and read the result from the display.

9.3.1 Emptying the Measuring Cell

1. Place the filling tube in a vessel which is suitable for sample disposal.
2. Empty the measuring cell by pressing the pump lever.

10 Cleaning and Storing the Instrument

10.1 Cleaning the Measuring Cell

Clean the measuring cell with a suitable solvent regularly before and after each measurement series to ensure the long-term accuracy of your results. If the measuring cell is not sufficiently cleaned, residue may form in the measuring cell and lead to inaccurate measurement results. Depending on the application, also cleaning between measurements may be required.

NOTICE

Make sure the solvent you use for cleaning is suitable. The Snap 40 pump cylinder is not resistant to aggressive cleaning liquids such as acetone or methyl ethyl ketone (MEK).

NOTICE

Do not use any mechanical action for cleaning the measuring cell.

1. Empty the measuring cell (see Chapter 9.3.1).
 2. Fill the measuring cell with a suitable solvent.
 3. Pump the solvent through the whole measuring system several times.
 4. Empty the measuring cell.
-

NOTICE

When only cleaning with one cleaning liquid, do not use tap water instead of distilled water. Otherwise, scale in the measuring cell may falsify your measured result.

10.2 Cleaning the Filling Pump

Clean the filling pump regularly, to ensure a smooth gliding of the pump piston within the pump cylinder.

1. Open the pump lock and remove the pump.
2. Rinse the pump under running tap water while moving the pump piston up and down.
3. Dry the pump with a lint-free cloth.
4. Insert the pump and close the pump lock.
5. Empty the measuring cell (see Chapter 9.3.1).

10.3 Cleaning the Pump Housing

Snap 40 is completely sealed against the outside and the pump. Therefore, you can rinse the pump housing with running tap water if the housing is dirty.

1. Open the pump lock and remove the pump.
2. Hold the pump housing under running water and rinse it thoroughly.
3. Insert the pump and close the pump lock.
4. Empty the measuring cell (see Chapter 9.3.1).

10.4 Cleaning the Housing and Display

To clean the housing and display, use a soft cloth dipped in ethanol or warm water. If necessary, you can use a mild solvent ($\text{pH} < 10$).

10.5 Storing the Instrument

Before storing the instrument for a longer period of time, clean the measuring cell as described in Chapter 10.1. This prevents any liquid residues freezing or drying in the measuring cell, which could destroy the cell or lead to measuring inaccuracies.

For storage lasting less than one day, fill the measuring cell with deionized water or solvent. If you have filled the liquid via the plastic syringe, leave the syringe in the adapter to stop the liquid running out.

11 Maintenance

11.1 Changing the Batteries

You can see the status of the battery by looking at the battery symbol in the top right-hand corner of the display.

The following status are possible:

Battery is full	
Battery is a little low	
Battery is low	
Battery is empty	

When the batteries are almost empty, the warning "Low Battery" is displayed. After the second warning "Battery Empty", the instrument switches itself off.

If the batteries are flat, you can easily replace them.

WARNING

When opening the battery compartment or exchanging batteries, sparks may be generated which can cause an explosion or fire in hazardous areas. Serious injuries are possible.

- Never open the battery cover in hazardous areas.
- Only replace the batteries outside hazardous areas.

NOTICE

- Always observe the correct polarity when replacing the batteries.
- Only use batteries of the same type and with the same level of battery charge.

1. Empty the measuring cell (see Chapter 9.3.1).
2. Open the pump lock and remove the pump.

Fig. 11 - 1 Removing the pump

3. Screw out the screw of the pump lock using an Allen key size 2.5 and remove the screw.
4. Remove the pump lock by sliding it outwards.

Fig. 11 - 2 Removing the pump lock

5. Use the Allen key size 2.5 to open the battery cover.

Fig. 11 - 3 Removing the battery cover

6. Remove the flat batteries by turning the instrument upside down.
7. Insert the new batteries. The picture on the back of the instrument shows the correct positioning of the batteries.
8. Close the battery cover using the Allen key.
9. Insert the pump lock again and slide it back to the original position. Make sure you slide the pump lock in the correct direction.
10. Fix the pump lock again using the Allen screw.
11. Insert the pump again.
12. Close the pump lock.

The instrument is ready for operation.

13. Check that the instrument is working properly by switching it on and checking the display.

11.2 Carrying out a Firmware Update

You can import the new firmware version into the instrument via the integrated IrDA interface. You need a PC with an infrared interface or a PC with an IrDA USB adapter installed and connected. Contact your Anton Paar representative if you need an IrDA USB adapter for your PC.

To establish a connection to a PC

1. If your PC has no IrDA interface, install the IrDA USB adapter on your PC.
2. Hold the IrDA interface of Snap 40 to the IrDA interface of your PC to establish the connection between the two instruments.

To import the firmware

1. Establish a connection between the Snap 40 and the PC as described above.
2. Press <Menu> and select **"Setup > Firmware Update"**.
3. On the PC, right mouse-click on the file you want to send and select **"Send to > A nearby computer"**.

After the successful file import, Snap 40 checks the validity of the update file.

4. Press <OK> to start the installation when "Import valid" is shown.

The progress bar shows you the progress of the installation. After successful installation the instrument switches itself back on.

11.3 Device Information

Information about the instrument is saved in the instrument and can be accessed at any time. This information includes the following:

- Manufacturer
 - Instrument name
 - Serial number
 - Firmware version
 - Date of the firmware version
 - Bootloader version
 - Date of the bootloader version
 - Hardware status
- Press <Menu> and select **"Information"**.

The device information is displayed.

12 Readjustment

If the measuring cell is not sufficiently cleaned, residue can form in the measuring cell. This can be noticed if the density values for deionized water begin to deviate. Therefore, perform regularly check measurements with deionized water.

If the check measurement shows a deviation of more than $\pm 0.04\%$, a readjustment is required. A readjustment is only possible if the calculated difference between the measured and the theoretical value is smaller than $2\% \text{ v/v}$ and the temperature of the water is between 15 and $25\text{ }^{\circ}\text{C}$ (59 and $77\text{ }^{\circ}\text{F}$).

12.1 Check Measurement

1. Clean the measuring cell as described in Chapter 10.1 and rinse the measuring cell with deionized water until no more solvent residues are present.
2. Fill the measuring cell with deionized water at a temperature between 15 and $25\text{ }^{\circ}\text{C}$ (59 and $77\text{ }^{\circ}\text{F}$).
3. Press the  key or press <Menu> and select **"Adjustment > Water Check"**.
4. Confirm with <OK> that the measuring cell has been filled with deionized water. The deviation of the measured density value compared to the reference density value is shown in %.
5. If the check measurement shows a deviation higher than $\pm 0.04\%$, perform a readjustment.

12.2 Water Adjustment

One option is to readjust using deionized water.

1. Clean the measuring cell as described in Chapter 10.1 and rinse the measuring cell until no more solvent residues are present.
2. Press <Menu> and select **"Adjustment > Water Adjustment"**.
3. Fill the measuring cell with deionized water at approx. $20\text{ }^{\circ}\text{C}$ (15 to $25\text{ }^{\circ}\text{C}$) without bubbles.
4. Confirm with <OK> that the measuring cell is filled with deionized water.

The density value of the filled deionized water is now measured and compared with the adjustment data. The display shows the calculated deviation (Delta) in %.

5. Press <OK> to activate the new adjustment.

12.3 Restoring to Factory Adjustment

If you want to undo all your adjustments, you can restore the factory adjustment.

1. Press <Menu> and select "**Adjustment > Factory Adjustment**".

The display shows: "Are you sure that you want to reset to factory adjustment?"

2. Press <OK> to restore the factory adjustment.

The factory adjustment is restored.

13 Error Messages and Warnings

Error message	Cause
Out of Specification	The measured sample temperature is outside the specifications.
Out of Range	- The measured value is outside the specified range. Possible reasons: - The temperature is too high or too low. - The density is too high or too low. - The measured value is invalid. - The value of deionized water, measured during readjustment, is outside the allowed limits (2 % w/v).
Temperature Range	The temperature of the reference liquid for readjustment is outside the temperature specifications (+15 to +25 °C).
Criterion	The density of deionized water at 20 °C measured during readjustment deviates from factory adjustment by more than 2 % w/v.
Low Battery	The battery is almost empty.
Battery Empty	The battery is empty.
Error 01	The temperature is invalid.
Error 02	The period is invalid.
Error 03	Snap 40 has no adjustment data.
Error 04	Calculation error.

NOTICE

If the error messages "Error 01" to "Error 04" occur, send in your instrument for service.

14 Technical Data

Measuring range:	
Alcohol:	0 to 100 % v/v
Temperature:	5 to 30 °C (34.1 to 86 °F)
Accuracy:	
Alcohol:	0.2 % v/v
Temperature:	0.2 °C (0.4 °F)
Repeatability, s. d.:	
Alcohol:	0.1 % v/v
Temperature:	0.1 °C (0.2 °F)
Resolution:	
Alcohol:	0.01 % v/v
Temperature:	0.1 °C (0.1 °F)
Sample volume:	approx. 2 mL
Sample temperature:	0 to 100 °C (32 to 212 °F)
Ambient temperature^a:	0 to +35 °C (32 to 95 °F)
Storage temperature^b:	-20 to +70 °C (-4 to 158 °F)
Humidity:	5 to 90 % relative air humidity, non-condensing
Protection class:	IP54 (use in light rain or snow possible)
Display:	LCD monochrome, graphical, 41 x 55 mm, 160 x 240 pixel
Data memory:	no data storage
Interfaces:	InOBEX/IrLPT
Power supply:	
Battery type:	Alkaline battery 1.5V, type EN91 (LR06, AA) from Energizer Industrial
Battery life:	> 100 hours
Dimensions^c:	140 x 138 x 27 mm (5.5 x 5.4 x 1.1 inches)
Weight (incl. batteries):	345 g (12.2 ounces)

a The sample must not freeze in the oscillator.

b Empty the measuring cell before storing the instrument.

c Without pump and filling tube

15 Wetted Parts

The following materials are in contact with samples and the cleaning liquids.

Material	Part
Borosilicate glass	Measuring cell, pump cylinder
PP (Polypropylen) Faradex MS002	Housing
PTFE (Polytetrafluoroethylene)	Filling tube, piston for hand pump
PVDF (Polyvinylidene fluoride)	Connection block, bolt
Viton	O-ring

Appendix A: Density of Water¹ (0 °C to 40 °C)

T °C	0.0	0.1	0.2	0.3	0.4	0.5	0.6	0.7	0.8	0.9
0	.99984	.99985	.99985	.99986	.99987	.99987	.99988	.99988	.99989	.99989
1	.99990	.99990	.99991	.99991	.99992	.99992	.99993	.99993	.99993	.99994
2	.99994	.99994	.99995	.99995	.99995	.99995	.99996	.99996	.99996	.99996
3	.99996	.99997	.99997	.99997	.99997	.99997	.99997	.99997	.99997	.99997
4	.99997	.99997	.99997	.99997	.99997	.99997	.99997	.99997	.99997	.99997
5	.99996	.99996	.99996	.99996	.99996	.99996	.99995	.99995	.99995	.99994
6	.99994	.99994	.99993	.99993	.99993	.99992	.99992	.99991	.99991	.99991
7	.99990	.99990	.99989	.99989	.99988	.99988	.99987	.99987	.99986	.99985
8	.99985	.99984	.99984	.99983	.99982	.99982	.99981	.99980	.99980	.99979
9	.99978	.99977	.99977	.99976	.99975	.99974	.99973	.99973	.99972	.99971
10	.99970	.99969	.99968	.99967	.99966	.99965	.99964	.99963	.99962	.99961
11	.99960	.99959	.99958	.99957	.99956	.99955	.99954	.99953	.99952	.99951
12	.99950	.99949	.99947	.99946	.99945	.99944	.99943	.99941	.99940	.99939
13	.99938	.99936	.99935	.99934	.99933	.99931	.99930	.99929	.99927	.99926
14	.99924	.99923	.99922	.99920	.99919	.99917	.99916	.99914	.99913	.99911
15	.99910	.99908	.99907	.99905	.99904	.99902	.99901	.99899	.99897	.99896
16	.99894	.99893	.99891	.99889	.99888	.99886	.99884	.99883	.99881	.99879
17	.99877	.99876	.99874	.99872	.99870	.99869	.99867	.99865	.99863	.99861
18	.99859	.99858	.99856	.99854	.99852	.99850	.99848	.99846	.99844	.99842
19	.99840	.99838	.99836	.99835	.99833	.99831	.99828	.99826	.99824	.99822
20	.99820	.99818	.99816	.99814	.99812	.99810	.99808	.99806	.99803	.99801
21	.99799	.99797	.99795	.99793	.99790	.99788	.99786	.99784	.99781	.99779
22	.99777	.99775	.99772	.99770	.99768	.99765	.99763	.99761	.99758	.99756
23	.99754	.99751	.99749	.99747	.99744	.99742	.99739	.99737	.99734	.99732
24	.99730	.99727	.99725	.99722	.99720	.99717	.99715	.99712	.99709	.99707
25	.99704	.99702	.99699	.99697	.99694	.99691	.99689	.99686	.99683	.99681
26	.99678	.99676	.99673	.99670	.99667	.99665	.99662	.99659	.99657	.99654
27	.99651	.99648	.99646	.99643	.99640	.99637	.99634	.99632	.99629	.99626
28	.99623	.99620	.99617	.99615	.99612	.99609	.99606	.99603	.99600	.99597
29	.99594	.99591	.99588	.99585	.99582	.99579	.99577	.99574	.99571	.99568
30	.99564	.99561	.99558	.99555	.99552	.99549	.99546	.99543	.99540	.99537
31	.99534	.99531	.99528	.99524	.99521	.99518	.99515	.99512	.99509	.99506
32	.99502	.99499	.99496	.99493	.99490	.99486	.99483	.99480	.99477	.99473
33	.99470	.99467	.99463	.99460	.99457	.99454	.99450	.99447	.99444	.99440
34	.99437	.99433	.99430	.99427	.99423	.99420	.99417	.99413	.99410	.99406
35	.99403	.99399	.99396	.99393	.99389	.99386	.99382	.99379	.99375	.99372
36	.99368	.99365	.99361	.99358	.99354	.99350	.99347	.99343	.99340	.99336
37	.99333	.99329	.99325	.99322	.99318	.99314	.99311	.99307	.99304	.99300
38	.99296	.99292	.99289	.99285	.99281	.99278	.99274	.99270	.99267	.99263
39	.99259	.99255	.99252	.99248	.99244	.99240	.99236	.99233	.99229	.99225
40	.99221	.99217	.99214	.99210	.99206	.99202	.99198	.99194	.99190	.99186

1 Excerpt from: Spieweck, F. & Bettin, H.: Review: Solid and liquid density determination. Technisches Messen 59 (1992), pp. 285–292

Appendix B: Measuring Unit Details

ID	Measuring unit	Minimum	Maximum	Unit	Format
1	Alcohol % v/v @ 20 °C	000.00	100.00	% Vol	nnn.nn
2	Alcohol US @ 60 °F	000.00	200.00	Proof	nnn.nn

ID	Measuring unit	P_{min} [g/cm ³]	P_{max} [g/cm ³]	T_{min} [°C]	T_{max} [°C]	T_{min} valid [°C]	T_{max} valid [°C]
1	Alcohol % v/v @ 20 °C	0.5	1.0	-10	40	5	30
2	Alcohol US @ 60 °F						

Appendix C: CE Declarations of Conformity

CE Declaration of Conformity

Anton Paar

Anton Paar GmbH hereby declares that the product listed below in the version offered for sale meets all the basic requirements of the applicable sections of the relevant EU directives in design and type.

This declaration will be deemed invalid should any unauthorized modifications be made to the product. Follow the information given in the instruction manual when setting up and operating the instrument.

Product designation: **SNAP 40 Alcoholmeter**

Model: **SNAP 40**

Manufacturer: **Anton Paar GmbH**

The product meets the requirements of the following directives:

- **Electromagnetic Compatibility 2004/108/EC**

Applied standards:

EN 61325-1:2006

Electrical equipment for measurement, control and laboratory use
- EMC requirements - Part 1: General requirements

The product is classified as a class B equipment and is intended for the use in industrial area.

- **Low Voltage Directive 2006/95/EC**

Applied standards:

EN 61010-1:2010

Safety requirements for electrical equipment for measurement,
control and laboratory use Part 1: General requirements

Ing. Peter Kettisch
Division Manager

DI Günter Holer
Department Manager

Appendix D: Firmware Versions

Firmware version	Release date	Document number	Comments
V1.00	30.04.2014	C96IB005ML-A	First released version.

Appendix E: Menu Tree

Adjustment	Water Check	
	Water Adjustment	
	Factory Adjustment	
Unit	% Vol; Proof	
Setup	Display Contrast	-8, ..., +8
	Backlight	On; Auto; Off
	Energy Saving	Off; 3 min; 5 min; 10 min
	Beep	On; Off
	Language	English; German; French; Spanish
	Firmware Update	
Information		

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Herausgegeben von Anton Paar GmbH, Gedruckt in Österreich. Copyright © 2014
Hersteller: **Anton Paar GmbH**