
Messen

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

Automatic drift control AutoRead:

<AR>  Switch on the drift control.

<RUN/ENTER>  AR flashes. Wait until the measured value is stable and AR stops flashing. You can cancel AutoRead at any time using <RUN/ENTER>.


Salinity correction (concentration measurement only):

<▲> Switch on the salinity correction. Sal appears on the display.


<▼> Switch off the salinity correction. Sal disappears from the display.


Calibration

Calibration in water vapor saturated air, e. g. in an air calibration vessel.

<CAL>, until  Put the D. O. probe in the air calibration vessel. The sponge in the vessel must be moist (not wet!).

<RUN/ENTER>  AutoRead measurement.

 As soon as the measured value is stable, the relative slope appears.

 Return to the measuring mode using <M/O₂>.

Entering the salinity of the measuring solution:

<CAL>, until  Using <▲> <▼>, set the salinity. Range: 0.0 ... 70.0 mg/kg Return to the measuring mode using <M/O₂>.

Calibration via a comparison measurement.

<CAL>, until  Immerse the D. O. probe in the calibration solution.

<RUN/ENTER>  AutoRead measurement.

 As soon as the measured value is stable, the measured value and the correction factor for the relative slope appear.

Using <▲> <▼>, adjust the measured value to the value of the comparison measurement.


Press <M/O₂> to return to the measuring mode.


Display indication:

The unit of the measured value flashes when the adjusted correction factor is $\neq 1$.



Storing

Storing the displayed measured value manually:


<STO>  No.: Number of the next free storage location.

<RUN/ENTER>  Ident: ID number of the data record. Using <▲> <▼>, change the ID number. Finish storing with <RUN/ENTER>.

Automatic storing function AutoStore:

<RUN/ENTER> +  Using <▲> <▼>, change the storing interval. <STO>  (OFF = Storing function switched off).

<RUN/ENTER>  Number of the storing locations that are still free.

<RUN/ENTER>  Ident: ID number of the data record. Using <▲> <▼>, change the ID number. Start AutoStore using <RUN/ENTER>.

Data record: Storage location number, date/time, measured values, temperature, temperature measuring procedure, ID number.

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
**Note to this
operating manual**

This operating manual contains a description of

- all basic functions,
- all instructions for a safe operation, and
- all technical data in a compact form.

A more detailed description with notes for special applications is available as a pdf document via internet under <http://www.schott.com/labinstruments>.

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handylab pH 11 - List of contents

Technical data	14
Safety	16
Display and jack field	16
Switching on the measuring instrument	17
Measuring.	17
Calibration	18
Configuration	22
Resetting to default settings	22

Technical data

Dimensions and weight	Length [mm]	172
	Width [mm]	80
	Height [mm]	37
	Weight [kg]	Approx. 0.3

Mechanical structure	Type of protection	IP 66
-----------------------------	--------------------	-------

Electrical safety	Protective class	III
--------------------------	------------------	-----


Test certificates	cETLus, CE
--------------------------	------------

Ambient conditions	Storage	- 25 °C ... + 65 °C
	Operation	-10 °C ... + 55 °C
	Relative atmospheric humidity	< 90 % of annual average

pH/ORP measuring ranges	Measuring range		Resolution
	pH	-2.000 ... +19.999 - 2.00 ... + 19.99	0.001 0.01
U [mV]	- 999.9 ... + 999.9 - 1999 ... + 1999	0.1 1	
T [°C]	- 5.0 ... + 105.0	0.1	

Precision of pH/ORP (± 1 digit)	pH (± 2 pH units from the calibration point)	± 0.005 at + 15 °C ... + 35 °C ± 0.01
	U [mV]	± 0.3 at + 15 °C ... + 35 °C ± 1

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
FCC Class A Equipment Statement

Note: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

Safety

Safety instructions

The individual chapters of this operating manual use safety instructions such as the label shown below to indicate various hazards or dangers:


Caution


indicates instructions that must be followed precisely in order to avoid the possibility of slight injuries or damage to the instrument or the environment.

Display and jack field

Display

Jack field

- 1 pH electrode
- 2 Reference electrode
- 3 Temperature probe

Caution

Only connect electrodes to the measuring instrument that do not return any unallowed voltages or currents (> SELV and > current circuit with current limiting). Almost all electrodes - in particular SCHOTT electrodes - fulfill these conditions.

Indication	Asymétrie [mV]	Pente [mV/pH]
Calibration incorrecte	< -30 ou bien > 30	< -62 ou bien > -50

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
Switching on the measuring instrument

- Press the $\langle \text{ON} \rangle$ key.
The test screen appears briefly on the display.
After this, the measuring instrument automatically switches to the measuring mode. The display shows the relevant measured value.

Measuring

Overview of the measuring modes:

Special functions:

Selecting the pH resolution

If you press $\langle \text{pH} \rangle$ while pressing the $\langle \text{RUN/ENTER} \rangle$ key, the resolution switches between 0.01 and 0.001.

AutoRead (drift control)

The AutoRead function checks the stability of the measurement signal (except for the measurement of the ORP voltage). Activate AutoRead with $\langle \text{AR} \rangle$. The current measured value is frozen. Thus, the $\langle \text{AR} \rangle$ key can also be used as a "hold function" to freeze the measured values. Press $\langle \text{RUN/ENTER} \rangle$ to start the AutoRead measurement. During the AutoRead measurement, AR flashes on the display until a stable measured value is reached. This can be terminated at any time and the current measured value taken over with $\langle \text{RUN/ENTER} \rangle$.

Temperature determination for pH measurements

You can use pH electrodes with and without temperature sensors. The measuring instrument recognizes which temperature sensors are connected and switches automatically to the correct mode for the temperature measurement (display shows TP).
Manual temperature input: Set the temperature value with $\langle \blacktriangle \rangle$ $\langle \blacktriangledown \rangle$.

Calibration

AutoCal TEC

This function is specially adjusted to SCHOTT technical buffer solutions as a fully automatic **one, two, or three point calibration** procedure. The buffer solutions are automatically recognized by the measuring instrument. Depending on the instrument setting, the instrument displays the relevant buffer nominal value or the current electrode voltage in mV.

Valid buffer (values at 25 °C): 2.00 / 4.01 / 7.00 / 10.01

Note

The calibration at pH 10.01 is specially optimized for the SCHOTT technical buffer solution, pH 10.01. Other buffer solutions can lead to an erroneous calibration. The correct buffer solutions are given in the SCHOTT catalog or in the Internet.

AutoCal DIN

This function is like AutoCal TEC, but matched to buffer solutions according to DIN 19266.


Valid buffer (values at 25 °C): 1.679 / 4.006 / 6.865 / 9.180

ConCal

This function is a conventional **two-point calibration** using two buffer solutions (pH 7.0 ± 0.5 and any other buffer solution) or a **single-point calibration** using any buffer solution and is used as a high-speed method.

Calibration evaluation

After the calibration, the measuring instrument automatically evaluates the calibration. The asymmetry and slope are evaluated separately. The worst case evaluation appears on the display.


Display	Asymmetry [mV]	Slope [mV/pH]
	-15 ... +15	-60.5 ... -58
	-20 ... +20	-58 ... -57
	-25 ... +25	-61 ... -60.5 or -57 ... -56
	-30 ... +30	-62 ... -61 or -56 ... -50
Clean the electrode according to the electrode operating manual		
 Invalid calibration	< -30 or > 30	< -62 or > -50

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Configuration

The following setting can be adapted individually to your requirements (the delivery status is shown in bold typeface):

Display during pH calibration

Buffer nominal value, or current electrode voltage

Procedure for configuring (factory settings appear in bold typeface):

- Switch off the measuring instrument.
- Press the <⏻> key while pressing the <pH> key.
- Select the required display during the pH calibration with <▲> <▼>. Selection: **Buffer nominal value (pH)** or current electrode voltage (mV).
- Then press the <RUN/ENTER> key. The measuring instrument automatically switches to the measuring mode.

Resetting to default settings

Default settings

The following settings are reset to the status on delivery when a reset is performed (initialized):

pH/ORP (pH In)

Measuring mode	pH
Slope	-59.16 mV/pH
Asymmetry	0 mV
Manual temperature input	25 °C
Resolution (pH display)	0.01
Display during pH calibration	Buffer nominal value

Reset procedure:

- Press <CAL> while pressing the <RUN/ENTER> key.
- Press <RUN/ENTER> to reset the measurement parameters. The instrument changes to the measuring mode.

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Copyright